

InStride Shoes was founded to provide comfortable, well fitting shoes for walking and leisure activities that provide significant relief for people with foot problems ranging from bunions, hammertoes and corns, to more significant diabetes related problems.

Developed in conjunction with Podiatrists, Pedorthists, and other medical and shoe construction specialists, InStride shoes provide both a protective function as well as relief from regular footwear due to their ability to provide a more precise and comfortable fit.



Seamless Linings

Seamless Dri-Lex® linings with Aegis® anti-microbial shield provide a dry foot environment and comfortable fit.



Dri-Lex®

Dri-Lex® fabrics are designed to manage moisture within a shoe environment to keep the foot surface dry and comfortable. To move moisture, all Dri-Lex® fabrics combine fiber content, yarn features and fabric construction with hydrophilic (water loving) components.



Pre-molded Heel Counter

Pre-molded counter for added stability in the heel area.




Shank Support

Shank support for stability. Airport friendly.



Tylite Gel Pad

Tylite gel heel pad provides cushioning when walking and standing.



Stride Guide™

Stride Guide™ bisects outsole and midsole to provide lateral and medial stability as well as shank support.



Technology



InStride Shoes feature the **CFS** – Custom Fit System – made with a series of removable insoles that provide the extra depth that allow InStride to fit even the most difficult feet with comfort beyond compare. InStride shoes are made with high quality full-grain leather, with Nubuc and Lycra styles available with lace or strap closures. When combined with non-irritating seamless linings, slip-resistant outsoles, the latest in shock absorbing and cushioning technology, InStride shoes assure the best possible combination of quality construction and good looking, contemporary styling for everyone.

InStride Custom Fit System

The InStride Custom Fit System (CFS) provides 2–3 levels of comfort and meets the changing needs of the “hard to fit feet”. The 2–3 removable insoles can be used in any combination to customize the shoe or allow sufficient room for an orthotic insert. Meets Medicare TSD guidelines. Accommodates new “K”-code inserts – A5512.



Full Grain Leather

Genuine high quality full grain leather uppers.



Stretch Lycra Vamp

For added comfort (not shown on display shoe)



Stridelite Midsole

STRIDELITE midsole for lightweight comfort or compression molded EVA midsole.



Toe Box

Generous rounded reinforced toe for all day comfort.



Trampoline Action innersole reciprocates energy expended on each step and promotes forward motion.

Sure Flex grooves in outsole allow for optimal flexibility.

Slip Resistant Outsoles

InDurance outsole provides high abrasion/high mileage durability or Dupont Everlon® slip resistant compound in outsole provides traction and skid resistance on slippery surfaces.

