

Drew Last Recommendations for Foot Problems:

Women's Styles:

Pedic 381

- *Deep oblique toe shape provides more forefoot capacity
- *8/8 heel height



Pedic 18

- *Snug fitting level heel seat
- *Sculpted high arch
- *Short back last



Pedic 88

- *A universal last used in many applications for extreme foot conditions
- *Deep, broad, oblique, forepart offers more capacity than any of our other women's lasts
- *Wide shank area provides superior lateral stability – extremely versatile for orthotic applications
- *8/8 heel height



Pedic 16

- *Added depth combined with broad forepart
- *Wide shank area for superior lateral stability
- *Excellent for applications to treat minor foot conditions
- *6/8 heel height



Pedic 46

*Wide shank area provides superior lateral stability and also makes this last extremely versatile for orthotic applications

*Deep, broad, oblique, forepart offers maximum capacity

*6/8 heel height offers excellent balance and weight distribution



Pedic 200

*Snug fitting heel seat

*Short back last

*Excellent for applications to treat minor foot conditions

*6/8 heel height



Pedic 246

*Snug fitting heel with displacement in cuboid

*Wide shank area provides excellent lateral stability

*Excellent last to accommodate prescribed orthotics

*8/8 heel height

*Additional overall capacity



Pedic 48

*A universal last used in many applications for extreme foot conditions

*Deep, broad, oblique forepart offers maximum capacity

*Short back last

*8/8 heel height



Men's Styles

Pedic 66

- *Snug fitting heel with extra displacement in cuboid
- *Wide shank area provides excellent lateral stability
- *Excellent last to accommodate prescribed orthotics
 - *6/8 heel height
- *Additional overall capacity



Pedic 466:

- *Snug fitting heel with displacement in cuboid
- *Wide shank area provides excellent lateral stability
- *Excellent last to accommodate prescribed orthotics
 - *6/8 heel height
- *Additional overall capacity

