

# **ComfortFit's New Pre-fabricated Orthotic Program; Designed For Cost Effective Treatment for *All Walks Of Life....***

*With today's strained economy, more and more patients are finding custom made orthotics financially out of reach. To combat this problem ComfortFit, along with its quality business partners, has developed a new and exciting line of affordable prefabricated orthotics. You can inventory these new products or just buy them as required for your patients.*



**I. Pre-fabricated orthotics for general use with athletic, casual and most dress shoes. Anti-Microbial, Anti-Fungal, High Wicking, Very Comfortable, Low Friction and Durable**



**ControlTecc - Boxed**

This new product features a firm VEVA molded shell with excellent rear and midfoot control. Special arch “peg” design provides control while minimizing weight and bulk. These new orthotics fit into most footwear and are ideal for heavy use and sports applications.

This new and innovative product comes complete with low friction and comfortable BambooTecc top covers that is an Ecco friendly fiber that has anti-odor, anti-microbial, anti-fungal and high wicking properties.



**SoftShell – Boxed**

Manufactured with an extremely soft VEVA molded shell with additional support, these light weight insoles have moderate control and a special shock absorbing gel heel plug. Perfect for the average patient as a replacement insole. This device will fit in almost all footwear styles.

This new insole is manufactured with a SilverTecc top cover that uses a new and exclusive technology from Samsung Corporation. Silvix provides a superior combination of high silver ionic content and fiber. This new Nano Silver/Fabric is anti-microbial and anti-fungal while also enhancing the physical properties of comfort. Soft surface with low friction.

## II. Functional Pre-fabricated Orthotics for Men's and Women's Sport and Dress shoes



**Comfort Dress**

A functional, intrinsically posted, full width, 3/4 length semi-rigid foot orthotic covered in a sulcus length soft tissue top cover for easy fit and comfortable feel. This is a product of choice when a generic, cost effective solution is desired by your patient in a shoe with limited room.



**Comfort Active**

An orthotic for maximum control, extrinsically posted and covered in 1/8" of soft tissue top covering. With firmer shell and a deeper heel cup than our Dress model, this insole creates the perfect cost effective sport device.

### III. Pre-fabricated Orthotics for Children



An exclusive ComfortFit product that is a cost effective alternative for children who are constantly outgrowing their custom molded orthotics. This device provides highly effective therapy for children needing control of severe pronation and ankle instability.

A patented multi-layer, firm anti-pronatory orthotic that will control the subtalar joint. Firm medial and lateral flanges with a foam Pelite soft top cover and anti-slip silicon bottom.

### IV. Pre-fabricated orthotics for the diabetic patient or those with very sensitive feet



A prefabricated orthotic designed for patients with diabetes, arthritis, or any other conditions that result in sensitive feet. These Plastazote covered insoles also provide superb resilience, cushioning and shock absorption. Medicare approved: suggested billing code A0628



Heat moldable insoles that provide maximum protection for the diabetic and sensitive foot. Full length multi-density insole with a Plastazote top cover designed to be used with all diabetic shoe types. Medicare approved: suggested billing code A5512