

INTRODUCING

# THE *Richie* GAUNTLET™

*Legitimate Varus/Valgus control of the hindfoot*



**RICHIE GAUNTLET WITH  
LATERAL ARCH SUSPENDER**



**RICHIE GAUNTLET WITH  
MEDIAL ARCH SUSPENDER**

## FEATURES:

Varus or Valgus lifting strap • Intrinsic balanced AFO shell • Extrinsic rearfoot post

## INDICATIONS:

Rigid, non-reducible Adult Acquired Flatfoot (Stage III & IV) • Severe DJD or Deformity of hindfoot • Charcot Arthropathy

## CASTING INSTRUCTIONS:

Cast for the Richie Gauntlet using the STS Mid Leg Sock



Non-Weight & Weight bearing casts are accepted

Available in 2 colors, Tan & Chocolate



2 shell heights: 7" (left) and 9" (right)

## BILLING CODES:

L1940 AFO, Custom-fabricated  
L2275 Varus/Valgus Control  
L2280 Molded inner boot  
L2820 Soft interface

## AVAILABLE THROUGH:

**ComfortFit Orthotic Labs**  
**Roselle, New Jersey**  
**1-888-523-1600**  
**1-888-362-2021**